

*On your lenses?
In your lungs!*

Dent AirVac
ORAL AEROSOL EVACUATION EQUIPMENT

Contact Us Today!
877-336-8247
dentairvac@yahoo.com

DAV VII Turbo + Oral Aerosol Evacuation System

“At Source”

Only \$1095

- For Mercury Vapor, Ozone, BPA dust, abrasion powder, laser smoke, VOCs, airborne viruses...
- Disposable intake liners for moisture control
- Internal stiffener arm for extended non contact hose support
- 5 foot long, 4" - 3" - 5.5" diameter hands free intake hose with 5 inch laminar flow opening
- 4 stage filtration for moisture, particulate, odor and vapor control
- Compact, only 12" X 12" X 15", lightweight, only 35 pounds
- Rugged powder coated scratch resistant metal housing
- Powerful yet quiet 770 cfm, 70 dbl, 110v, 60hz vacuum
- Casters for multiple operatory mobility
- Prefilter for particulate capture down to 10 micron
- HEPA filter for particulate capture down to 0.3 micron
- Granulated activated carbon filter for odors and vapors
- Variable Speed Motor Controller Option
- Limited Lifetime Warranty

